

	Environmental Analysis Teaching and Research Laboratory	Date: X/XX/XXXX	Number: X
	Standard Operating Procedure	Title: Accumet XL600 Meter	
	Approved By: TBD	Revision Date: November 29, 2017	

1. Scope and Application

1.1 This SOP describes how to use the Accumet meter system to measure and monitor a variety of water-based electrochemical parameters including pH, conductivity, temperature, salinity, dissolved oxygen, and ion-selective mV.

2. Summary of Method

2.1 This SOP does this ...

Contents

1	Scope and Application	1
2	Summary of Method	1
3	Definitions	3
4	Biases and Interferences	3
5	Health and Safety	3
	Safety and Personnel Protective Equipment	3
6	Personnel & Training Responsibilities	3
7	Required Materials and Apparati	3
8	Reagents and Standards	4
9	Estimated Time	4
10	Sample Collection, Preservation, and Storage	4
11	Procedure	4
	pH	4
	Temperature	4
	BOD (Biological Oxygen Demand)	4
	Conductivity	4

12 Data Analysis and Calculations	5
Graphing	5
13 QC/QA Criteria	5
14 Trouble Shooting	5
No response or all buffers read the same pH	5
Slow electrode response with excessive crystallization inside the electrode	5
Slow electrode response with due to clogged junction	6
Dried Salt Deposits	6
Slow electrode response, noisy, unstable or erratic readings	6
15 References	6

3. Definitions

3.1 Term 1: Electrode- The glass electrodes are ion-selective meaning the glass membrane is sensitive to a specific ion. For example, the pH electrode is sensitive to hydrogen ions.

3.2 Term 2: Reference Cavity-

3.3 Term 3: Buffer- A buffer solution is one which resists changes in pH when small quantities of an acid or an alkali are added to it. Buffer solutions are used as a means of keeping pH at a nearly constant value in a wide variety of chemical applications.

4. Biases and Interferences

4.1 Biases and interferences can come from ...

5. Health and Safety

5.1 Describe the risk ...

Safety and Personnel Protective Equipment

6. Personnel & Training Responsibilities

6.1 Researchers training is required before this the procedures in this method can be used...

6.2 Researchers using this SOP should be trained for the following SOPs

- SOP01 Laboratory Safety
- SOP02 Field Safety

7. Required Materials and Apparati

7.1 Accumet XL600 Benchtop Meter (XL94105590)

7.2 BOD Probe (2446146 445)

7.3 Temperature Probe (93X306506 505)

7.4 pH Electrode (13 620 183A)

7.5 Conductivity Electrode (13 620 100)

7.6 Standards for calibration (pH 4.00, 7.00, 10.00, conductivity, and a small bottle that can be used to make the DO standard)

8. Reagents and Standards

9. Estimated Time

9.1 This procedure requires XX minutes ...

10. Sample Collection, Preservation, and Storage

11. Procedure

pH

11.1 Calibration

11.2 The level of electrolyte in the outer cavity should be kept above the level of the solution being measured to prevent reverse electrolyte flow. The electrolyte need only be immersed far enough to cover both the glass pH sensing bulb and reference junction to obtain accurate readings.

11.3 Electrode Diagram

Temperature

11.4 Calibration

BOD (Biological Oxygen Demand)

11.5 Calibration

Conductivity

11.6 Calibration

The accumet electrodes utilize a simple twist open and close design. The purple electrodes are Tris compatible and are less prone to cloggin by sulfides, heavy metals or proeints. A Blue band indicates a single junction electrode designed for general use in camples without sulfides, heavy metals, or proteins.

12. Data Analysis and Calculations

Graphing

- 12.1** Real-time data can be viewed on a graph to view changes over brief or extended periods. Time is plotted in seconds. The graph refreshes every hour from the start of graphing.
- 12.2** Select Show Graph. A graph will be displayed on the lower half of the screen where various measurement details were previously displayed.
- 12.3** To track live measurement data on the graph, select Start Plotting. The measurement data will continue until Stop Plotting is selected. When stopped, the graph can be dragged left/right and up/down.
- 12.4** Select the Zoom In or Zoom Out options to examine the data more closely or broadly.
- 12.5** To remove the graph display and revert to display measurement data on the lower half of the screen, select Hide Graph. Note: Measurement data will continue to be tracked on the graph until Stop Plotting is selected.

13. QC/QA Criteria

14. Trouble Shooting

No response or all buffers read the same pH

- 14.1** Verify that the correct meter input and channels are connected. Try removing the electrode storage bottle or rubber bulb guard to ensure contact with sample. If the meter has automatically frozen reading, be sure the meter's hold or auto read feature are off. Lastly, try replacing the electrode.

Slow electrode response with excessive crystallization inside the electrode

- 14.2** Verify that the fill hole is open during measurements and closed during storage. Otherwise, the electrolyte flow could be clogged with supersaturated electrolyte. Flush and refill electrode. Remove the filling solution through the fill hole with a syringe or by shaking it upside down. Repeatedly flush and rinse the reference cavity with clean water at a temp of 60-80 degrees celsius to dissolve crystals. Replace the filling solution and apply gentle pressure to filling hole. Hydrate the electrode with in storage solution of pH 4 buffer.

Slow electrode response with due to clogged junction

- 14.3** For protein deposits, prepare a 1% pepsin solution in 0.1M HCl and soak the reference junction for one hour. Rinse the electrode with distilled water. For general cleaning, heat a diluted KCl solution to 60 to 80 degrees Celsius and soak the sensing bulb for about 10 minutes. Allow the electrode to cool in unheated KCl solution.

Dried Salt Deposits

- 14.4** Dissolve the salt deposits in warm tap water and then soak the electrode briefly in pH 4 buffer.

Slow electrode response, noisy, unstable or erratic readings

- 14.5** Clean the electrode with mild detergent and warm water and then hydrate. Allow the electrode to reach sample temp. Take 30 minute readings and soak the electrode in pH buffer for 1 minute between measurements.

15. References

- 15.1** APHA, AWWA, WEF. (2012) Standard Methods for examination of water and wastewater. 22nd American Public Health Association (Eds.). Washington. 1360 pp. (2014).

- 15.2** Online Instruction Manual can be found here: [link](#)