

	Environmental Analysis Teaching and Research Laboratory	Date: 02/02/2017	Number: 20
	Standard Operating Procedure	Title: Ozone Generator	
	Approved By: Marc Los Huertos	Revision Date: February 2, 2017	

## 1. Scope and Application

**1.1** The scope of this SOP is train researchers in how to effectively use the Laboratory Ozone Generator

**1.2** As a researcher, occasions will arise in which fabricated ozone gas must be used in order to test a hypothesis. Through using this generator, researchers can create ozone gas in small quantities, from dry air or oxygen and with negative or positive pressures.

## 2. Summary of Method

## 3. Definitions

## 4. Operation

**4.1** In order to operate the ozone generator there are several procedures that must be followed.

- 1. Electrically connect the unit to the mains supply utilising the cable supply lead supplied with the unit.
- 2. Start feeding gas through the generator and set the generator to the required flowrate.
- 3. Depress the main ON-OFF switch on the generator front panel which will illuminate indicating that ozone gas is being produced.
- 4. Set the variable control knob to the output required by utilising the output graphs included in this manual.
- 5. The unit takes ten minutes initially to reach its normal operating temperature and output.
- NOTE: The red FAULT indicator will illuminate briefly after pressing the main ON-OFF switch due to the fractional delay in the power board relays latching.

## **5. Interferences**

## **6. Health and Safety**

### **6.1. Safety and Personnel Protective Equipment**

## **7. Personnel & Training Responsibilities**

Researchers training to use the Eosense chambers and Picarro analyzer include the following components:

Researchers using this SOP should be trained for the following SOPs:

- SOP03 Field Work
- SOP04 Electrical Power in the Field

## **8. Required Materials**

### **8.1. Item 1 w/catalog number!**

### **8.2. Item 2**

## **9. Estimated Time**

**9.1** This will take XX minutes...

## **10. Procedure**

**10.1** Prepare ...

**10.2**

## **11. References**

**11.1** APHA, AWWA, WEF. (2012) Standard Methods for examination of water and wastewater. 22nd American Public Health Association (Eds.). Washington. 1360 pp. (2014).